



Technology®
ABclonal

Leader in Biomolecular
Solutions for Life Science

ABclonal Adjustable Volume Pipette

Backed by an excellent R&D center, state-of-the-art manufacturing facilities certified by ISO 9001 & ISO 13485, and an end-to-end quality management system, we deliver superior quality products that offer unrivaled price-performance ratio.

Product List

Cat.No.	Description	Size
AI20010P	0.5-10μL Single-channel Adjustable Volume Pipette	1Set
AI20200P	20-200μL Single-channel Adjustable Volume Pipette	1Set
AI21000P	100-1000μL Single-channel Adjustable Volume Pipette	1Set



Product Features

- 1 Ergonomic design of low plunger pressure avoids repetitive strain injuries;
- 2 3 volumes cover a range of 0.5μL to 1000μL;
- 3 Large 4-position volume display and high display accuracy up to 0.01μL;
- 4 Durable high-tech synthetic material provides excellent resistance to high temperature, corrosion and weather;
- 5 Reasonable structure and exquisite technology make flexible assembling and disassembling;
- 6 Tip ejector makes easy and quick tip change;
- 7 High-elastic tip cone provides good seal and easy removal of tips;
- 8 Easy to calibrate and maintain with tools supplied;
- 9 Calibrated in accordance with ISO8655 and each pipette supplied with individual test certificate;
- 10 The lower section can be autoclaved.

DY-C-22JUL07-EN01